



AN 385

Analysis of Contaminants by DHS/GC/MS

May 7, 2007 (Version 3.0)

Discussion

Dynamic Head Space/Gas Chromatography/Mass Spectrometry combines the advantages of utilizing three techniques for solving difficult problems in organic analysis. Head Space Analysis (HSA) was initiated for analysis of volatile compounds in matrices that could not be directly injected into a gas chromatograph. Dynamic Headspace is a non-equilibrium technique.

The technique works by subjecting the sample to elevated temperatures for a period of time to drive volatile compounds from the sample matrix into the atmosphere above the sample, called the "head space". The inert gas helium continuously sweeps the chamber, carrying with it the components that have outgassed. The purge gas is passed through an adsorbent trap that contains (in order): glass beads, tenax™, and charcoal and is cooled down to -100°C. At the conclusion of the experiment, the trap is thermally desorbed (in reverse flow) directly into the inlet of a gas chromatograph. The components that were present on the trap are separated by Gas Chromatography and identified using Mass Spectrometry. An external standard usually deuterated n-hexadecane is used to provide semi-quantitative results.

This short note presents an analysis of good disk drives and failed disk drives. The breather filters and internal absorbers were analyzed by Dynamic Headspace GC/MS to identify contaminants that might be causing the failures. Figure 1 compares a siloxane (retention times 13.66, 16.16) contaminated breather filter from a failed drive with a non-siloxane contaminated breather filter from a good drive. Figure 2 compares a siloxane (retention times 11.01, 13.66 and 16.16) contaminated inner absorber from the same failed drive with a non-siloxane contaminated inner absorber from the same good drive.




Figure 1. Comparison of Breather filters

DHS/GC/MS analysis of this type is applicable to any type of sample, not just the disk drive filters used in this example. Any sample that fits in the DHS chambers (1 1/2" diameter x 4" long) can be outgassed at any temperature between 45°C and 300°C, normally for a time period of one hour to three hours. Any volatile organic components will be detected and semi-quantitative results will be provided. Detection limits are as low as 10 ng/component.




Figure 2. Comparison of Inner Absorbers

United States Locations

Tempe, Arizona
 +1 480 239 0602 info.az@eaglabs.com
 +1 602 470 2655 fax

Sunnyvale, California
 810 Kifer Road
 +1 408 530 3500 info.ca@eaglabs.com
 +1 408 530 3501 fax

1135 E Arques Avenue
 +1 408 738 3033
 +1 408 738 3035 fax

785 Lucerne Drive
 +1 408 737 3892
 +1 408 737 3916 fax

Peabody, Massachusetts
 +1 978 278 9500 info.ma@eaglabs.com
 +1 978 278 9501 fax

Chanhassen, Minnesota
 +1 952 828 6411 info.mn@eaglabs.com
 +1 952 828 6449 fax

East Windsor, New Jersey
 +1 609 371 4800 info.nj@eaglabs.com
 +1 609 371 5666 fax

Syracuse, New York
 +1 315 431 9900 info.ny@eaglabs.com
 +1 315 431 9800 fax

Raleigh, North Carolina
 +1 919 829 7041 info.nc@eaglabs.com
 +1 919 829 5518 fax

Round Rock, Texas
 +1 512 671 9500 info.tx@eaglabs.com
 +1 512 671 9501 fax

International Locations

Shanghai, China
 + 86 21 6879 6088 info.cn@eaglabs.com
 + 86 21 6879 9086 fax

Tournefeuille, France
 + 33 5 61 73 15 29 info.fr@eaglabs.com
 + 33 5 61 73 15 67 fax

Frankfurt, Germany
 + 49 (0) 693053213 info.de@eaglabs.com
 + 49 (0) 69307941 fax

Tokyo, Japan
 + 81 3 5396 0531 info.jp@eaglabs.com
 + 81 3 5396 1930 fax

HsinChu, Taiwan
 + 886 3 5632303 info.tw@eaglabs.com
 + 886 3 5632306 fax

Uxbridge, United Kingdom
 + 44 (0) 1895 811194 info.uk@eaglabs.com
 + 44 (0) 1895 810350 fax