




AN 426

SIMS Detection Limits of Selected Elements in ZnO Under Normal Depth Profiling Conditions


May 7, 2007 (Version 2.0)

O ₂ /SIMS Positive Secondary Ion Detection		Cs/SIMS Negative Secondary Ion Detection	
Element	DL (atoms/cm ³)	Element	DL (atoms/cm ³)
B	1E+14	H	2E+17
Mg	5E+13	C	1E+16
Al	5E+13	N	2E+16
K	1E+13	F	5E+15
Ti	1E+14	Si	1E+17
Cr	5E+14	P	2E+15
Fe	5E+14	As	1E+16
Ni	1E+15		
Cu	1E+15		
Ga	5E+14		
In	1E+15		
Cd	5E+16		

Example Profiles Provided on Back Side



Overlay of two separate SIMS profiles of an arsenic implant in ZnO. Excellent reproducibility is demonstrated in this example as well as the count-rate-limited detection limits of approximately $1\text{e}16$ at/cm³.



Overlay of multiple species in a single SIMS profile for a ZnO sample containing implants of Ti, Cr, Ni and Cu.

United States Locations

Tempe, Arizona
+1 480 239 0602 info.az@eaglabs.com
+1 602 470 2655 fax

Sunnyvale, California
810 Kifer Road
+1 408 530 3500 info.ca@eaglabs.com
+1 408 530 3501 fax

1135 E Arques Avenue
+1 408 738 3033
+1 408 738 3035 fax

785 Lucerne Drive
+1 408 737 3892
+1 408 737 3916 fax

Peabody, Massachusetts
+1 978 278 9500 info.ma@eaglabs.com
+1 978 278 9501 fax

Chanhassen, Minnesota
+1 952 828 6411 info.mn@eaglabs.com
+1 952 828 6449 fax

East Windsor, New Jersey
+1 609 371 4800 info.nj@eaglabs.com
+1 609 371 5666 fax

Syracuse, New York
+1 315 431 9900 info.ny@eaglabs.com
+1 315 431 9800 fax

Raleigh, North Carolina
+1 919 829 7041 info.nc@eaglabs.com
+1 919 829 5518 fax

Round Rock, Texas
+1 512 671 9500 info.tx@eaglabs.com
+1 512 671 9501 fax

International Locations

Shanghai, China
+ 86 21 6879 6088 info.cn@eaglabs.com
+ 86 21 6879 9086 fax

Tournefeuille, France
+ 33 5 61 73 15 29 info.fr@eaglabs.com
+ 33 5 61 73 15 67 fax

Frankfurt, Germany
+ 49 (0) 693053213 info.de@eaglabs.com
+ 49 (0) 69307941 fax

Tokyo, Japan
+ 81 3 5396 0531 info.jp@eaglabs.com
+ 81 3 5396 1930 fax

HsinChu, Taiwan
+ 886 3 5632303 info.tw@eaglabs.com
+ 886 3 5632306 fax

Uxbridge, United Kingdom
+ 44 (0) 1895 811194 info.uk@eaglabs.com
+ 44 (0) 1895 810350 fax