

# Analytical Solutions

for BioTechnology

May 9, 2007 (Version 2.0)

## BN 1422


### Evaluation of Cleaning Efficacy


#### Discussion

The removal of surface contamination is a common goal of all cleaning processes. Unfortunately all cleaning steps are not equally effective, and some processes can actually leave behind additional contaminants. Electron Spectroscopy for Chemical Analysis (ESCA) is a surface analysis tool that is well qualified to characterize cleaning processes. ESCA's ability to quantify the concentration of all elemental species except hydrogen and helium makes it possible to monitor the efficiency of contamination removal as a function of cleaning processes.

The spectra shown below are from a glass sample before and after several cleaning steps. Peaks for sodium (Na) and carbon (C) decrease in intensity after cleaning, while peaks from the glass substrate (e.g. silicon and oxygen) increase in relative intensity as the surface gets cleaner. The table shows the concentration of the elements before and after cleaning. The first step reduces the level of Na, but increases the C concentration. The subsequent process step removes a majority of the C, leaving a cleaner glass surface.

	Atomic Concentrations (%)							
	Si	O	Na	C	N	Ca	Mg	F
As-Received	13	42	10	29	2.9	1.7	0.6	0.7
Clean I	19	44	0.2	32	3.7	0.5	0.5	0
Clean II	24	62	0.6	12	0.6	1.1	0.6	0





---

#### United States Locations

Tempe, Arizona  
+1 480 239 0602 [info.az@eaglabs.com](mailto:info.az@eaglabs.com)  
+1 602 470 2655 fax

Sunnyvale, California  
810 Kifer Road  
+1 408 530 3500 [info.ca@eaglabs.com](mailto:info.ca@eaglabs.com)  
+1 408 530 3501 fax

1135 E Arques Avenue  
+1 408 738 3033  
+1 408 738 3035 fax

785 Lucerne Drive  
+1 408 737 3892  
+1 408 737 3916 fax

Peabody, Massachusetts  
+1 978 278 9500 [info.ma@eaglabs.com](mailto:info.ma@eaglabs.com)  
+1 978 278 9501 fax

Chanhassen, Minnesota  
+1 952 828 6411 [info.mn@eaglabs.com](mailto:info.mn@eaglabs.com)  
+1 952 828 6449 fax

East Windsor, New Jersey  
+1 609 371 4800 [info.nj@eaglabs.com](mailto:info.nj@eaglabs.com)  
+1 609 371 5666 fax

Syracuse, New York  
+1 315 431 9900 [info.ny@eaglabs.com](mailto:info.ny@eaglabs.com)  
+1 315 431 9800 fax

Raleigh, North Carolina  
+1 919 829 7041 [info.nc@eaglabs.com](mailto:info.nc@eaglabs.com)  
+1 919 829 5518 fax

Round Rock, Texas  
+1 512 671 9500 [info.tx@eaglabs.com](mailto:info.tx@eaglabs.com)  
+1 512 671 9501 fax

#### International Locations

Shanghai, China  
+ 86 21 6879 6088 [info.cn@eaglabs.com](mailto:info.cn@eaglabs.com)  
+ 86 21 6879 9086 fax

Tournefeuille, France  
+ 33 5 61 73 15 29 [info.fr@eaglabs.com](mailto:info.fr@eaglabs.com)  
+ 33 5 61 73 15 67 fax

Frankfurt, Germany  
+ 49 (0) 693053213 [info.de@eaglabs.com](mailto:info.de@eaglabs.com)  
+ 49 (0) 69307941 fax

Tokyo, Japan  
+ 81 3 5396 0531 [info.jp@eaglabs.com](mailto:info.jp@eaglabs.com)  
+ 81 3 5396 1930 fax

HsinChu, Taiwan  
+ 886 3 5632303 [info.tw@eaglabs.com](mailto:info.tw@eaglabs.com)  
+ 886 3 5632306 fax

Uxbridge, United Kingdom  
+ 44 (0) 1895 811194 [info.uk@eaglabs.com](mailto:info.uk@eaglabs.com)  
+ 44 (0) 1895 810350 fax